

Size Measurement and Molecular Weight Estimation of Globins



Introduction

The oligomeric state of globins is a crucial piece of information to understand their structure-function relationships. It is known to vary significantly within the superfamily: proteins consisting of 1, 2, 4, 12, 24 globin subunits have been found in nature. Within the group of globins, it is particularly critical to assess the size of hemoglobin-based blood substitutes, as this has been shown to correlate with the extent of both extravasation and renal filtration, two factors responsible for the toxicity of these products once they are administered to patients.

The assessment of the oligomeric state of globins is usually carried out by gel filtration chromatography, a fairly time consuming technique. Moreover, chromatography measurements imply dilution of the sample which can complicate the data interpretation as concentration-dependent dissociation of the subunits can occur.

Dynamic light scattering (DLS) is a rapid batch technique allowing the measurement of samples without the need for dilution, hence avoiding the problem of possible complex dissociation. The technique is non-invasive and the samples can be recovered after measurements for further use.

DLS measures the time dependent fluctuations in scattered light from molecules in a solution to determine the translational diffusion coefficient and subsequently determine the hydrodynamic radius (R_H) from the Stokes-Einstein equation. The given R_H corresponds to the radius of a sphere diffusing with the same speed as the molecule in solution. The R_H

depends thus on the shape of the molecule, with an increased R_H for elongated molecules compared to spherical molecules of the same molecular weight. For similarly shaped molecules, a calibration curve can be constructed to estimate the molecular weight based on the measured R_H .

All globins share common fold and very similar 3D structures. In this application note, a calibration curve for globins has been constructed by measuring the following globins, where the molecular weight and subunit conformation is known: Myoglobin, Human hemoglobin, AHb1, AHb2 (two globins from *Arabidopsis thaliana*), rice Hb1 and Hemoglobin Polytaur (an auto-polymerizing recombinant globin proposed as a blood substitute). To extend the range of sizes measured,

two other proteins were added: Ferritin from spinach and Hemocyanin from *Eurypelma californicum*, as they are high molecular weight proteins with a globular shape. Table 1 summarizes the molecular weight and number of subunits in each globin.

The resulting calibration curve was thereafter used to estimate the apparent molecular weight of two hemoglobin-based oxygen carriers (HBOCs). Both of them are PEGylated hemoglobins, one called PEGHb^{oxy} and mimicking the MP4 product (Hemospan®) from Sangart, and the other is called PEGHb^{deoxy} [2]. Both are hemoglobin molecules (64 kDa) conjugated with a number of 6-8 polyethylene glycol (PEG) chains (5.5 kDa).





Figure 1: Intensity size distributions for Myoglobin, Rice Hb1, Human hemoglobin (HbA), recombinant hemoglobin (Polytaur) and hemocyanin



Methods

All measurements were performed on the Zetasizer Nano ZS, a dynamic light scattering (DLS) instrument, using non-invasive backscatter (NIBS) detection.

The calibration curve was constructed by plotting the measured hydrodynamic radius versus the known molecular weights. The data points were fitted with an equation of the following form $R_H = (a \cdot M_w)^b$ [1] to obtain a relationship between the hydrodynamic radius and the molecular weight (M_w) for globins.

All protein concentrations were in the range 5-500 μ M in saline buffer, where, for each protein, no subunit dissociation is predicted to take place.

All samples were filtered before measurement using 0.22 μ m or smaller pore size to remove large aggregates.

Results

Polydisperse Samples

Most of the samples showed a multimodal correlation function indicating the presence of several populations in the samples.

Oligomeric State of Globin Samples

Although the samples are polydisperse, the size distribution analysis provides information about the individual populations and allows for the determination of the oligomeric state of the globin samples. Figure 1 shows the intensity size distributions for a number of the globin samples. The figure clearly shows that the different samples have multiple size peaks.

The fraction of the sample in each size peak was obtained from the volume size distributions, determined using Mie theory. It was shown that, for all globin samples, the smallest size peak was the main component of

Table 1: Summary of the number of subunits, known molecular weight and peak statistics for various proteins

Protein	No. of Subunits	Known Mw [kDa]	Peak Mean [nm]	%Peak Intensity	%Peak Volume
Myoglobin	1	17.6	2.403	97.7	100
AHb2	1	17.8	2.341	53.2	99.9
Rice Hb1	2	37	2.985	90.9	100
AHb1	2	36	3.016	80.8	100
HbA	4	64.5	3.305	100	100
Polytaur	12	186	6.778	88.3	99.8
Ferritin	-	440	8.74	82.8	99.9
Hemocyanin	-	1720	12.36	100	100

the sample, making up at least 99.8% of the sample volume (Table 1).

The size of the main component in each sample was therefore determined as the mean value of the smallest size peak. In Table 1, a summary of the peak statistics for the main peak of each sample is shown. It is noteworthy that proteins of similar

molecular weight, such as Myoglobin (17.6kDa) and AHb2 (17.8kDa) have very similar hydrodynamic radii, 2.40nm and 2.34nm respectively. The same can be seen for Rice Hb1 (37kDa) and AHb1 (36kDa), where the R_H was measured to be 2.99nm and 3.02nm, respectively.





Figure 2: Hydrodynamic radius versus molecular weight for the various proteins measured (red circles). The experimental values have been fitted with a power function to obtain a calibration curve to estimate the molecular weight of unknown globins



These data confirm that it is possible to use DLS as a quick and user friendly screening tool to determine the oligomeric state of a globin sample, as well as a tool to detect the presence of aggregates in samples.

Hydrodynamic Size Versus Molecular Weight

The hydrodynamic radius of the main component in each globin sample is plotted versus the known molecular weight of the samples in figure 2. The data fit, giving the relationship between molecular weight (M_w) and hydrodynamic radius (R_H) has been added to figure 2.

The hydrodynamic radius data from hemocyanin from *Eurypelma californicum*, has been excluded from the fit, as the large measured size can be affected by angular dependence of the measurement.

The M_w of PEGylated Hemoglobins

The PEGylated hemoglobins contain multiple populations, with the smaller size peak being the main component of both samples. Figure 3 shows the intensity and volume size distributions for the two samples. The two samples differ as the PEGHb^{deoxy} sample shows a narrow main size peak, with R_H of 5.85nm, while the PEGHb^{oxy} sample contains at least two populations in its main size peak as indicated by a shoulder, with R_H values corresponding to 4.89nm and ca 9nm respectively.

Using the calibration curve previously established, the estimated apparent molecular weight of the main components of PEGHb^{oxy} are 108 kDa and 450kDa while the single peak of the PEGHb^{deoxy} sample corresponds to 165kDa. As the PEG chains conjugated to the protein affect the diffusion of the molecules, the calibration curve should only be used to estimate the corresponding weight of a globin that diffuses with the same speed as the PEGylated hemoglobin





Figure 3: (A) Intensity size distributions for PEGHboxy (red) and PEGHbdeoxy (green) (B) Volume size distributions for the sample samples. These measurements were made at 25°C.



products, hence the wording “apparent molecular weight”. The importance of measuring the size and thermal stability of these two products is discussed in more detail in another application note [2].

Conclusions

The size distributions of various globins have been determined using DLS. DLS allows the determination of the hydrodynamic radius of globin samples, and the identification of monomeric, dimeric and higher oligomers of globin samples. It has been shown that DLS can be effectively used to estimate the Mw and determine the oligomeric state and the fraction of aggregates present in globin samples.

By establishing a calibration curve for similarly shaped proteins, DLS was used to estimate the molecular weight of globins from their hydrodynamic radius, thus constituting a fast alternative or complement to gel filtration chromatography.

References

[1] Claes P., Dunford M., Kenney A., and Vardy P. “ An on-line dynamic light scattering instrument for macromolecular characterization” in “Laser Light scattering in biochemistry” SE Harding, DB Sattelle, and VA Bloomfield (eds) pp 66-76, RSC (1992).

[2] “Size, oligomeric state and thermal stability of PEGylated hemoglobins”, MRK1294-01

Acknowledgements

These results were obtained by Drs. Stefano Bruno, Serena Faggiano, and Luca Ronda at University of Parma, Italy with the help of Alice Brun, Alfatest srl, Italy.

Malvern Instruments Ltd

Enigma Business Park • Grovewood Road • Malvern • Worcestershire • UK • WR14 1XZ
Tel: +44 (0)1684 892456 • Fax: +44 (0)1684 892789

Malvern Instruments Worldwide

Sales and service centers in over 50 countries for details visit www.malvern.com/contact

© Malvern Instruments Ltd 2009

more information at www.malvern.com