


# MASTERSIZER3000

smarter particle sizing



## PERFORMANCE




0.01µm to 3500µm  
(10nm to 3.5mm)



Dry Dispersion



Wet Dispersion



Accuracy better than 1.0%\*†



Reproducibility better than 1.0% variation\*



<10sec

\* Sample and sample preparation dependant

† Accuracy defined for the recovery of the mean size of a narrow log-normal distribution

## SOFTWARE



SOP comparison tool for consistent results between users



One-click configurable data support



Advanced measurement management functions



Analysis support



Upgradable

## GLOBAL SUPPORT



+44 (0) 1684 892 456



Engineering Support



Application Support



[salesinfo@malvern.com](mailto:salesinfo@malvern.com)



Webinars



eLearning Courses



[www.malvern.com/ms3000](http://www.malvern.com/ms3000)



1,000s of Citations